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Soil Sample (Rendzina Bv-Horizont, 150 mg)

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
Application Note

AN-S-270510-E-01

Determination of the CS content of cement using the rapid CS cube

Task

The analysis of cement samples is a difficult task for elemental analyzers due to the high calcium content in cement. Alkaline ions produced during the combustion process damage the quartz reaction tubes of the analyzer and disturb the S measurements. The rapid CS cube is able to accurately determine the CS content of cement using WO₃ and borax as additive.

Instrument	Sample
Basis: rapid CS cube	Quantity: 20 mg
Mode: CS	Consistency: solid
Periphery: solid sample former	Preparation: not necessary

Specification

All samples have been weighed into tin boats. A mixture of WO₃ powder and borax (2:1) has been added to the samples in order to bind earth alkaline ions and to prevent the formation of CaSO₄, which would bias the sulphur measurements. The WO₃/borax additive should be mixed well with the sample, and has to be added in a ratio 1:10 (sample:WO₃/borax).

Procedure

The instrument has been calibrated using certified soil standards.

Sample	C [%]	S [%]
Portland cement standard-1 (S = 1.85%)	0.23 ± 0.015	1.82 ± 0.036
Portland cement standard-2 (S = 0.82%)	0.25 ± 0.035	0.89 ± 0.071
Portland cement	0.67 ± 0.008	1.54 ± 0.036
Blast furnace cement	0.33 ± 0.028	2.28 ± 0.043
Portland blast furnace (slag) cement	0.61 ± 0.020	1.12 ± 0.028
hinder	5.90 ± 0.042	0.57 ± 0.028
Portland fly ash cement	1.42 ± 0.026	1.05 ± 0.054

Results

The results show that the sulphur contents of the two NIST cement standards have been measured with a very high accuracy. Furthermore, a very good precision has been found for all samples.

Due to the high combustion temperature of the rapid CS cube, WO₃/borax can be used as additive to cement and concrete samples, which efficiently increases the lifetime of the quartz tubes and leads to reliable sulphur measurements.

The rapid CS cube is very suitable for applications in the cement industry.

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