

INSTRUCTION FOR USE

PRODUCT GROUP OF ORTHODONTIC WAX



TECHNICAL DATA

Manufacturer:	BİLKİM TIBBİ ÜRÜNLER SAN.TİC.LTD.ŞTİ 1212/1 Sk. No:1/1 Menemen, Izmir/TURKEY
Product/Product Type:	Orthodontic Wax
Product Form:	Solid paraffinic wax with or without flavour
Material Type:	Fully Refined Paraffinic Wax

PRODUCT DESCRIPTION

Orthodontic wax is specially formulated to relieve pain, irritation or discomfort which may sometimes be caused by orthodontic devices such as the brackets and wires of braces, ligatures or even partial dentures which have wire clasps.

INSTRUCTION FOR USE

- A little peice of the wax is taken fromthe closed plastic case.
- The wax is lined up to 5 sections. Only one section is removed from the body of the wax, and rubbed over the braces to remove the irrugarities wich lead to relief.
- The rest of the wax is put back in the case to keep for usage a next time.

PRECAUTIONS

- Do not mix with other products.
- Do not use this product for purposes other than those indicated above.
- If symptoms of irritation and sensitivity are observed while using this product, discontinue use and seek medical treatment.

STORAGE

- Protect from heat and direct sunlight.
- Store in a dry place at room temperature (5°C -35°C / 41°F-95°F).
- Do not use after the expiration date.

Explanation of the Markings on the Packaging



Symbol for Manufacturer Information



Symbol for Manufacturing Date



Symbol for Lot Number

INSTRUCTION FOR USE
PRODUCT GROUP OF ORTHODONTIC WAX



Symbol for Catalogue Number



Symbol for Temperature Limits



Symbol for Expiration Date



Symbol for Keep Away from Direct Sunlight



Symbol for Follow the Instructions for Use



Symbol for Non-sterile Product



Confirmation: The product complies with the applicable European Directives